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## EN Instructions for use alphatech® rotary instruments

### 1.1 Brief description of rotary instruments

alphatech® surgical rotary instruments and auxiliary instruments are part of the alphatech® implant system. The rotary instruments belonging to the alphatech® implant system are drills used to prepare an osteotomy for a dental implant procedure with special regard to the alphatech® implant system. The alphatech® drills are available in different dimensions and shapes for the individual preparation of the osteotomy. All instruments are available in surgical kits, allowing the clinician to easily organise, store and transport the components. The drills are arranged in order of increasing diameter, according to the recommended drilling sequence, which facilitates orientation during implant placement. Colour-coded fields indicate the corresponding diameter.

### 1.2 Purpose

alphatech® rotary instruments and auxiliary instruments are intended for planning, preparing the osteotomy before placing the implant or for inserting implants or placing prosthetic components of the alphatech® implant system.

### Specific intended use

Mucosal punch	The mucosal punch is used to remove the gingiva before the first bone drilling.
Rose borer	Surface smoothing of the alveolar ridge after exposure of the cortical bone.
Cortical drill	Opening of the cortical area.
2.0 mm drill and depth stops	Expansion of the initial canal to Ø 2 mm and intended implant length. Depth stops for intended implant length.
2.8 mm drill bit and depth stops	Widening of the drill channel to Ø 2.8 mm according to the specified cavity depth. Depth stops for intended implant length.
Universal drill	According to the implant diameter and bone density.
Finisher	According to the implant diameter. Always use the finisher to be able to insert the implant neck into the cavity without stress and to avoid pressure necrosis on the implant.
Trepan drill	Trephine burs are indicated for the preparation of the bone flap after apical osteotomy. The bur enables the preparation of bone cylinders for the purpose of grafting.
Slim Line® 2.0 mm Drill	The 2 mm Slim-Line bone drill is indicated exclusively for the alphatech® Slim-Line implant variants. The drill is indicated for widening the drill channel to the corresponding implant depth and diameter.
Slim Line® 3-sided drill	The 3-sided bone drill is exclusively indicated for the alphatech® Slim-Line implant variants. The drill is indicated for cortical perforation, filigree tip and narrow atrophic alveolar ridges.

### 1.3 Target patient group

The target group is patients with chewing function restrictions due to missing or damaged teeth, for whom implant placement is carried out as a result.

### 1.4 Materials

The rotating instruments are made of rust-resistant medical stainless steel (EN standard 1.4197, ASTM 420F).

### 1.5 Indication

The rotary instruments and auxiliary instruments are intended for use in dental procedures performed to restore partially edentulous or edentulous patients with implants or prosthetic components of the alphatech® implant system. These instruments are indicated for intraoral use by a dentist. The cutting instruments are indicated for drilling or milling in the mandibular or maxillary bone during implant surgery procedures and can be used to prepare the jawbone and oral soft tissues.

**IMPORTANT:** Please use the universal drills according to the bone quality:  
2-ring drills for bone density D3 and D4  
3-ring drills for bone density D1 and D2

### 1.6 Contraindication

Contraindicated are patients who are medically unsuitable for oral surgery or who are allergic or hypersensitive to medical stainless steel or any of its alloy components.

### 1.7 Side effects and interaction

The clinical outcome of dental treatment is influenced by several variables. The following possible residual risks and side effects are associated with this product. These may require additional treatment in the dental office:

- Fracture of the implant/loss of the prosthetic component, problems with biting/chewing/speaking
- Bone compression
- Bone damage
- Injury to adjacent teeth/antagonists
- Complaints, discomfort
- Rotational vertigo
- Hypersensitivity reactions/allergic reactions
- Gum injuries
- localised pain
- longer wound healing/healing than expected
- Loss of the implant
- other toxicity reactions
- poor aesthetic result
- Surgical removal of the implant
- Ingestion/aspiration of small parts during the procedure
- Perforation of the maxillary sinus
- Nerve damage that can lead to chronic pain
- Paresthesia
- Dysaesthesia
- Swelling
- local inflammation
- Bruise
- Resorption of the upper/lower ridge
- Systemic or local infection (including peri-implantitis, periodontitis, gingivitis, fistula).
- minor bleeding.

### 2.1 General handling

The rotary instruments are intended for multiple use. However, the user (doctor) decides on the possibility of repeated use and the frequency of use of the instruments on the basis of the respective case of application and the possible wear of the products at his own responsibility. In case of doubt, the instruments should always be sorted out early and replaced by new ones. Rotating instruments must be clamped as deeply as possible and must be brought to the intended speed before being applied to the object. The optimal speeds to be used for the various rotary instruments can be found in the surgical manual or, in individual cases, on the label. Any liability is excluded if the instructions are disregarded. Avoid jamming or levering instruments, as this may lead to an increased risk of breakage. Depending on the application, the use of protective goggles is recommended when using rotating instruments. When working on dry materials, work should be carried out with an extractor. Particularly when using hand tools, care should be taken to use them gently and sensitively.

**IMPORTANT:** Thermal damage to the bone should be avoided in any case by optimal cooling and adherence to the recommended speeds: Perform all drilling for implant bed preparation at a speed of 400 – maximum 1,000 rpm and sufficient external cooling with sterile, physiological, cooled saline solution.

**IMPORTANT:** Avoid tilting during drilling, otherwise there is a risk of crestal, uncontrolled widening of the implant bed.

### 2.2 Dosage form

All rotary instruments are packed non-sterile and must be sterilised before use, depending on the intended use.

### 2.3 Cleaning and disinfection

Before first use on the patient and immediately after each use, all used instruments must be disinfected, cleaned and sterilised. If possible, a mechanical process (disinfector) should be used for cleaning and disinfection. A manual process – also using an ultrasonic bath – should only be used if a mechanical process is not available, due to the significantly lower effectiveness and reproducibility. Pre-treatment is mandatory in both cases.

### Pre-treatment

Pre-treatment must be carried out immediately after application (within max. 2 hours). During pre-treatment, coarse impurities must be removed from the products. Running water or a disinfectant solution should be used for this purpose. It must be ensured that the disinfectant is aldehyde-free, otherwise blood contamination may become fixed on the products. Furthermore, the disinfectant should have a tested efficacy (e.g. VAH/DGHH or FDA approval or CE marking), be suitable for disinfecting the products and be compatible with the products in terms of material resistance. Only use the nylon brushes or clean, soft and lint-free cloths intended for this purpose to remove contamination manually. Never use metal brushes or steel wool to remove contamination from the products! For products with cavities (e.g. trephine drills), rinse the cavities three times using a disposable syringe (minimum volume 5–10 ml) and a cannula. Please note that the disinfectant used in the pre-treatment is only for personal protection and cannot replace the disinfection step to be carried out later – after cleaning has been completed. Use by a dentist is indicated.

### A: Mechanical cleaning/disinfection (disinfector/cleaning and disinfection device)

1. Disassemble the product as much as possible.
2. In the case of products with cavities (lumens), ensure that all cavities have been effectively flushed and all coarse impurities have been removed as part of the pre-cleaning process.
3. Place the disassembled products in the disinfector. Make sure that the products do not touch each other. If possible: Connect all lumens of the products to the rinsing connection of the disinfector using a suitable rinsing adapter.
4. Start the programme.
5. Remove the products from the disinfector at the end of the programme.
6. Check and pack the cleaned and disinfected products as directly as possible after removal from the disinfector.

**IMPORTANT:** When selecting the disinfector, pay attention to the following:

- The disinfector must always have a tested efficacy (DGHH or FDA approval or CE marking according to DIN EN ISO 15883).
- Use of a tested programme for thermal disinfection (AO value > 3000 or at least 5 min. at 90°C for older units) – chemical disinfection is not to be used due to the risk of disinfectant residues on the products.
- Suitability of the programme used for the products with sufficient rinsing cycles
- Rinse with sterile or low-germ (max. 10 germs/ml) or low-endotoxin (max. 0.25 endotoxin units/ml) water.
- Drying process must be carried out with filtered air
- Regular maintenance and inspection of the washer-disinfector must be carried out.

**IMPORTANT:** When selecting the detergent system to be used, pay attention to the following:

- Cleaning agent system must basically be suitable for cleaning products made of metal and plastics
- If thermal disinfection is not used, a suitable disinfectant with tested effectiveness (e.g. VHA/DGHH or FDA approval or CE marking) must also be used, which is compatible with the cleaning agent used.
- The chemicals used must be fully compatible with the products in terms of material resistance
- The concentrations specified for cleaning agents and disinfectants must be strictly adhered to.

### B: Manual cleaning and disinfection

1. Dismantle the products as much as possible.
2. Place the disassembled products in the cleaning bath for the specified exposure time so that the products are sufficiently covered with liquid (if necessary, support with ultrasound treatment, but at least careful brushing with the nylon brushes provided). Make sure that the products do not touch each other. For products with cavities (lumen), rinse the cavities of the products five times at the beginning or at the end of the exposure time using a disposable syringe (minimum volume 5–10 ml) and a cannula.
3. Remove the products from the cleaning bath and rinse them thoroughly with water at least three times. For products with cavities (lumen), rinse all cavities five times at the beginning or end of the exposure time using a disposable syringe (minimum volume 5–10 ml) and a cannula.
4. Check the products.
5. Place the disassembled, cleaned and checked products in the disinfection bath for the specified exposure time so that the products are sufficiently covered with liquid. Make sure that the products do not touch each other. For products with cavities (lumens), rinse the cavities of the products five times at the beginning or end of the exposure time using a disposable syringe (minimum volume 5–10 ml) and a cannula.
6. Then remove the products from the disinfection bath and rinse them thoroughly with water at least three times. For products with cavities (lumens), rinse the cavities of the products five times at the beginning or end of the exposure time using a disposable syringe (minimum volume 5–10ml) and a cannula.
7. Dry the products by blowing them off with filtered compressed air.
8. Pack the products as soon as possible after removal.

**IMPORTANT:** When selecting the cleaning agents and disinfectants to be used, pay attention to the following:

- Basic suitability for cleaning or disinfection of products made of metals or plastics
- Suitability of the cleaning agent for ultrasonic cleaning – if applicable – without foam development
- Use of a disinfectant with tested efficacy (e.g. VHA/DGHH or FDA approval or CE marking), which is compatible with the cleaning agent used.
- The chemicals used must be fully compatible with the products in terms of material resistance.
- The concentrations of the detergents and disinfectants must be strictly adhered to.
- Combined cleaning agents and disinfectants should not be used if possible. Only in cases of very low contamination (no visible contamination) can combinable cleaning agents and disinfectants be used.
- Only freshly prepared solutions, only sterile or low-germ (max. 10 germs/ml) and low-endotoxin (max. 25 endotoxin units/ml) water may be used and only filtered air may be used for the drying process.

### 2.4 Control

After cleaning or cleaning/disinfection, visually inspect all products for corrosion, damaged surfaces/blank spots, dull or broken blades/chipping, damage to shape (e.g. bent, not round) and contamination and remove damaged products. Products that still show contamination and soiling despite cleaning/disinfection must be cleaned and disinfected again.

### 2.5 Reusability

The end of the product life of products for reuse is basically determined by wear and tear and damage from use. Therefore, it is necessary to check the products after or before each new use. If necessary, known restrictions on the frequency of use of products must be observed. If no wear is visible, the instruments can sometimes be used for longer. Please note the information on the frequency of use of the various instruments, as damaged and deformed blades can cause vibrations and lead to poor preparation edges and rough surfaces. Blunt and broken instruments result in high contact forces and thus increase the working temperature. This can lead to damage to the surrounding tissue (pulp). The guideline value for the frequency of use of the rotating instruments belonging to the alphatech® implant system, which may deviate from the actual service life depending on the application and / or material processed, is 4 times.

### 2.6 Maintenance

The disassembled products must be reassembled according to any specific instruction that may be available. Instrument oils must not be used!

### 2.7 Sterilisation

Only the sterilisation procedures listed below are recommended for sterilisation: Steam sterilisation:

- Fractionated vacuum process or gravity process (with sufficient product drying) – The gravity process is only permitted if the fractionated vacuum process is not available, as it is less effective.
- Steam steriliser according to DIN EN 13060 or DIN EN 285
- Validated according to DIN EN ISO 17665 (previously DIN EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ (picking) and product-specific performance assessment (PQ))
- Maximum sterilisation temperature 134°C or 273°F plus tolerance according to DIN EN ISO 17665 (previously DIN EN 554/ANSI AAMI ISO 11134)
- Sterilisation time (exposure time at sterilisation temperature) at least 20 minutes at 121°C or 250°F or at least 5 (or 18 minutes for prion inactivation) minutes at 132°C or 270°F/134°C or 273°F

### 2.8 Packing

Sort the cleaned and disinfected products into the associated drill stand or sterilisation tray/OP tray. Please pack the products or the drill racks/sterilisation trays/OP trays in single-use sterilisation packaging (single or double packaging) and/or sterilisation containers that meet the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilisation (temperature resistance up to at least 137°C or 279°F) with sufficient steam permeability
- Sufficient protection of the products or sterilisation packaging against mechanical damage.
- Regular maintenance according to the manufacturer's specifications (sterilisation container)
- For individual packaging: The packaging must be large enough so that the seal is not under tension.

### 2.9 Storage

Store the products in their original packaging in a dry place at room temperature and protected from direct sunlight. Further storage of the products should be in suitable, hygienically maintained containers (protected from moisture and recontamination), observing the maximum storage period. After sterilisation, the products must be stored dry and dust-free in the sterilisation packaging, whereby the expiry date of the sterile packaging resulting from the validation must be strictly observed. Improper storage can impair essential material properties and design characteristics and lead to functional failure.

### 2.10 Disposal

Disposal should be done in an environmentally sound manner and in accordance with local laws and regulations. Hazardous waste of contaminated products or sharps should be disposed of in appropriate containers that meet specific technical requirements.

### 3.1 Warnings

- During intraoral handling, all products must be secured against aspiration. Aspiration of products may result in infection or unforeseen injury to the patient.
- Do not use damaged, corroded or blunt instruments. Inspect the instruments before each use.
- During implant bed preparation and implant placement, the inferior alveolar nerve canal must be avoided. Nerve damage may result in anaesthesia, paraesthesia or dysaesthesia.
- To avoid bone necrosis, do not exceed the recommended insertion torques
- Ensure the correct drilling depth and use the depth marks on the drills, the depth stops and the alignment aids
- Use the drills and taps in clockwise rotation
- Use the drills in ascending diameters. Use an intermittent drilling technique and drill under ample cooling with pre-cooled sterile saline solution
- Bone quality should be taken into account when preparing the implant bed.

### 3.2 Safety instructions

- Unprotected touching of the instruments and parts by the user must be avoided at all costs (wear protective gloves).
- The instruments must only be used for the specified purpose. Always handle all instruments with extreme care to avoid needlestick injuries.
- Observe hygiene measures at the workplace according to the RKI guideline and carry out all work with powder- and latex-free gloves.
- Check the quality of the instruments regularly.
- When handling the components, be sure to follow the instructions described in the basic information on surgical and prosthetic procedures. Ensure sterile handling.
- Inspect the instruments before each use
- Always adhere to the recommended speeds.
- Use drilling templates and measuring probes, the depth marks on the drills, depth stops and alignment aids to ensure the correct depth and alignment of the hole.
- Insert the instruments into the bone without excessive force

Please observe the following instructions as part of your responsibility for the sterility of the products:

- Use of device- and product-specific validated procedures for cleaning, disinfection and sterilisation
- Regular maintenance and inspection of the equipment used (steriliser, disinfectant)
- Compliance with validated parameters for each cycle

When selecting cleaning agents and disinfectants, make sure that the following ingredients are not included:

- Organic, mineral and oxidising acids (minimum permissible pH value 5.5)
- Strong alkalis (maximum permissible pH value 8.5, neutral/enzymatic cleaner recommended)
- Organic solvents (e.g. alcohols, ethers, ketones, benzines), oxidising agents (e.g. hydrogen peroxides)
- Halogens (chlorine, iodine, bromine)
- Aromatic, halogenated hydrocarbons
- Never clean all products, drill stands and sterilisation trays/op trays with metal brushes or steel wool. All products and sterilisation trays/OP trays may only be exposed to temperatures up to a maximum of 137°C or 279°F!

### 4. Safety and liability

In case of non-observance of the instructions for use or warnings or in case of accidental misuse, the manufacturer accepts no liability for any damage that may occur. The reprocessor is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired results. This normally requires validation and routine monitoring of the procedure. Similarly, any deviation from the instructions provided should be carefully evaluated by the reprocessor for its effectiveness and possible adverse consequences. Serious incidents involving the instruments and implants should be reported to the national authorities as well as FMZ GmbH.

## 5. Explanation of symbols

	Manufacturer		Follow the instructions for use
	Date of manufacture		Attention, be careful when handling the appliance
	Usable until		Medical device
	Catalogue number		unique product identification
	Batch designation		Health Industry Bar Code (HIBC)
	Non-sterile		Prescription product
	Store dry		This product is a medical device according to MDR (EU) 2017/745 (Article 120).
	Protect from sunlight		

### Manufacturer

FMZ GmbH  
Charles-Darwin-Ring 3a I 18059 Rostock I Germany  
Phone +49 381 364 669-0  
E-mail: info@fmz-gmbh.de Internet: www.fmz-gmbh.de

Status 06/23, subject to change.  
GA-at-RI-300623-C