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**DESCRIPTION OF THE COMPONENTS**

The DIM analogs consist of a DIM analog body and a DIM knurled screw.

**DIM ANALOG COMPATIBILITY**

In each case, the DIM analog is compatible with the corresponding implant systems that are shown in Table 1.

Assignment of the DIM analogs to the implant systems:

Serie (DIM analogs)	Kompatibles Implantatsystem
ALP	FMZ alphatech® Tube-Line
BEG	BEGO Implant Systems/Semados®
BEGMU	BEGO Implant Systems Semados®MultiPlus - Abutment Niveau
BIO	BioHorizons® / internal
BPS	BEGO Implant Systems/Semados® SC / SCX / RS / RSX /RI-Implantate - Platform Switch Design
BRE	Bredent Medical® / Sky®
CAM	Altatec / Camlog®
COL	Altatec / Conelog®
E	Nobel Biocare™ / Replace Select®
EVMU	DENTSPLY Implants® / Astra Tech Implant System™ EV / Multi Unit
F	Nobel Biocare™ / Nobel Active™
H	Biomet 3i® / Certain®
HIO	HIOSSEN® / OSSTEM TS®
I	Biomet 3i® / Osseotite® External Hex®
ISY	Altatec / ISY®
K	Nobel Biocare™ / Brånemark®
L	Straumann / Bone Level®
MEG	MEGAGEN® / AnyRidge®
MEX	Medentis Medical® / ICX®
MIC	MIS® / C1®
MIS	MIS® / SEVEN®
N	Straumann / SynOcta®
NBMU	Nobel Biocare™ / Multi Unit
NEO	NEOSS Implant System® / NEOSS™
R	Zimmer Dental / Tapered Screw-Vent®
S	DENTSPLY Implants® / Astra Tech Implant System OsseoSpeed® TX
SEV	DENTSPLY Implants® / Astra Tech Implant System™ EV
SMU	DENTSPLY Implants® / Astra Tech Implant System OsseoSpeed® TX / Multi Unit
STRMU	Straumann® / Multi Unit
T	DENTSPLY Implants® / Xive®
TA	Tatum Implant System
THO	Thommen Medical® / SPI®
UPR	Sweden & Martina Implantology™ / Premium®
UTG	Sweden & Martina Implantology™ / Premium Kohno® / Premium TG®
Y	DENTSPLY Implants® ANKYLOS®
2-CON	Only in use with 2-CONnect Abutment BEGO Implant Systems/ Semados®; BioHorizons® / internal; Bredent Medical® / Sky®; Altatec® / Camlog®; Altatec / Conelog®; Nobel Biocare™ / Replace Select®; Nobel Biocare™ / Nobel Active™; Biomet 3i® / Certain®; HIOSSEN® / OSSTEM TS®; Biomet 3i® / Osseotite® External Hex®; Nobel Biocare™ / Brånemark®; Straumann / Bone Level®; MEGAGEN® / Medentis Medical® / ICX®; MIS® / C1®; MIS® / SEVEN®; Straumann / SynOcta®; NEOSS Implant System® / NEOSS™; Zimmer Dental / Tapered Screw-Vent®; DENTSPLY Implants® / Astra Tech Implant System OsseoSpeed® TX; DENTSPLY Implants® / Astra Tech Implant System™ EV; Sweden & Martina Implantology™ / Premium Kohno® / Premium TG®; DENTSPLY Implants® / Xive®; Thommen Medical® / SPI®

**MATERIALS**

DIM analog body: Ti6Al4V, Medical Grade 5, ASTM F136

DIM knurled screw: Stainless steel (1.4305).

**INDICATIONS**


DIM analogs are used to produce implant models in the digital workflow and also in the classic impression procedure using the DIM sleeve. DIM analogs make it possible to produce a prosthetic design and to check the fit.

## ACCESSORIES

- **DIM parallel countersink** (Stainless Steel 1.4305): By using a DIM parallel countersink, the opening to the knurled screw can be expanded and deburred when using DIM analogs; the positioning then takes place via the guide pin of this tool.
- **Sleeve** (Stainless Steel 1.4305): Through the placement of the sleeve on the analog body and fastening with the knurled screw, a rotation-free repositioning of the analog in the model is guaranteed.

## APPLICATIONS

### STEP-BY-STEP WORKFLOW: IN MODEL

Placement of the sleeve (2) on the analog body (1) and fastening with a knurled screw (3). The sleeve later remains in the model socket and ensures a rotation-free repositioning of the analog in the model.

Attachment of the DIM analog on the impression post in the impression by means of the impression post screw.

Fixing an extension sleeve (drinking straw, etc.) in place and covering the knurled screw with soft wax to prevent plaster from entering into the area between the sleeve and the screw.

When using socket systems/plates, cut the sleeves to length on plate contact and seal with wax.

Casting the impression and producing the gingival mask by using an impression made of modeling silicone.

The extension sleeves can now be removed from the model, the wax residue eliminated and, if necessary, the opening to the knurled screw for the holding fixture of the DIM parallel countersink guide pin can be exposed.

If necessary, the opening to the knurled screw can be expanded and deburred with the DIM parallel countersink; the positioning takes place via the guide pin of this tool.

By unscrewing the knurled screw, the analog can be removed; the sleeve remains in the model socket for the exact repositioning of the analog.

## STEP-BY-STEP WORKFLOW: DIGITAL

Start of the digital workflow: The design process is carried out and a suitable production process is controlled by means of CAM software to produce the model.

Positioning of the DIM analog in the model with an insertion tool: The three pins of the DIM analog are aligned to the placeholder of the model.

By using light pressure and gentle turning, the DIM analog is incorporated up to the perceptible stop in the model.

Positioned DIM analog: After removal of the insertion tool, the DIM analog is precisely anchored in the model. The knurled screw is fixed in place with the DIM analog.

Section view of the DIM analog fixed in the model: The nt-IQ Abutment library now makes it possible to precisely and reproducibly fabricate customized abutments in the digital workflow.

### APPLICATION INSTRUCTIONS

- The fixing in place of an extension sleeve and the covering up of the knurled screw must be carried out using soft wax to prevent plaster from entering into the area between the sleeve and the knurled screw and thus avoiding an inaccurate impression.
- When using a socket system/plate, the sleeve must first be cut to plate contact and sealed with wax in order to guarantee an accurate impression.
- Before use, check to ensure that the components are clean and free of damage.

### STORAGE

The components are to be stored dry and free of dust at room temperature.

nt-trading GmbH & Co KG / G.-Braun-Str. 18 / 76187 Karlsruhe, Germany / Tel: +49-721-91547160 /  
Fax: +49-721-91547161 / E-mail: info@nt-trading.com

Products indicated with ® are registered brand names of the manufacturers.

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Produkte gekennzeichnet mit © sind eingetragene Warenzeichen des entsprechenden Herstellers.