Hersteller:
FMZ GmbH
Charles-Darwin-Ring 3a
18059 Rostock, Germany

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Manufacturer:
FMZ GmbH
Charles-Darwin-Ring 3a
18059 Rostock, Germany

Manufacturer's instructions for use

Name of the medical device
BONITmatrix®

Composition and characteristics

BONITmatrix® is a synthetic, bioactive, completely absorbable material for the filling and reconstruction of bone defects. The production takes place in a patented procedure (DE 10003824). The material consists of up to 87 % by weight of two nano-crystalline calcium phosphates hydroxyl apatite (HA) and β -Tricalciumphosphate (β -TCP), which are embedded in a nano-porous silicon dioxide matrix (13 % by weight). The ratio from HA to β -TCP is 60/40. Since the material is synthetically manufactured, any transmission of infections or illnesses of animal origin is impossible. BONITmatrix® is characterized by an interconnecting pore system in the nano and micrometer range. The internal surface area is 90m²/g. This high porosity ensures a complete diffusion of body fluids (blood) and so favors the adsorption of Osteoinductive factors. When mixed high cohesive forces between the individual Granulate ensure that the augmentate remains in the defect. So-called micro motion is almost avoided and the regeneration of bone tissue can continue unimpaired. The material serves as a scaffold during bone formation.

Indications

- In Facio maxillary surgery, Periodontology and plastic surgery:
- for the filling and/or reconstruction of non-load bearing bone defects:
 - e.g.
 - Defects following removal of bone cysts,
 - Residual defects following osteotomy,
 - Other polymural alveolar defects,
 - Periodontal defects, also in combination with membranes,
 - Defects following extractions,
 - Defects after Apical Resection,
 - Extraction defects in preparation for implant therapy,
 - Sinus floor elevation,
 - Peri implant defects,
 - Defects after harvesting autologous bone.

Contraindications

BONITmatrix® may not be used with patients, who exhibit impeded wound healing or are receiving therapy which negatively influences the bone metabolism. In patients taking non-steroidal anti inflammatory agents (NSAIDS), as well as those with existing Gluco- and Mineral-corticosteroid therapy. A further contraindication is metabolic bone disorders.

Facio-maxillary surgery, plastic surgery, Periodontology:

- No specific contraindications apply. The general contraindications specified above apply.

Warnings

BONITmatrix® is not clinically tested on pregnant women, nursing mothers, children, at patients with metabolic bone diseases, with vascular disturbances, uncontrolled diabetes, Achondroplasia, anti-coagulation therapy, immunosuppressive therapy and therapy with Steroids. BONITmatrix® may not be implanted in acutely inflamed and infected regions. The area to be implanted must be free of soft and granular tissue. Contact with saliva or other liquids (e.g. physiological saline solution) must be avoided. Caution should be taken in patients with strong nicotine abuse.

Side effects

Implant related side effects e.g. wound healing disturbances, allergic reactions with proven allergy against individual implant components or inflammatory reactions have not been reported. Interactions with other medical products or medicines have not been reported in relation to BONITmatrix®. Concomitant medication: The use of analgesics is possible. The use of non-steroidal anti-inflammatory drugs (NSAIDS) can impede the regeneration of bone tissue.

Application

BONITmatrix® should be used exclusively by trained medical personnel with appropriate, professional training. The operational procedure depends upon the localization, type and size of the bone defect. For the success of BONITmatrix® a direct contact with vital bone tissue is absolutely necessary. Therefore the implantation of the bone regeneration material may take place only in vitalized defects. Prior to implantation it is necessary to mix BONITmatrix® into a malleable mass with the blood of the patient and if possible with autologous bone marrow. Autologous Spongiosa and/or bone chips may be in added to the BONITmatrix® mixture since they additionally support the bone new formation. In operations where local blood is difficult to obtain then venous blood should be used.

The implantation of BONITmatrix® into the defect should be carried out with a suitable instrument (e.g. Spatulas). BONITmatrix® mixture is form-stable and resistant to migration and therefore suitable for the treatment of larger defects. After the filling, closure of the bone defect is necessary, in order to prevent a washing out of the granules and the ingrowth of fibrous tissue. Ideally autologous Periosteum should be used. The augmented area should not be mechanically loaded during the initial healing phase. In joint-replacement procedures care should be taken to prevent BONITmatrix® from migrating into the articulating area. In load bearing indications internal or external stabilization techniques should be used where indicated, i.e. until BONITmatrix® is completely replaced by newly developed bone tissue. Wound drains should not be placed directly in contact with the bone substitute thereby preventing the evacuation of the cellular components so necessary for bone regeneration.

Dosage

The choice of size of the granulates depends upon the kind and size of the bone defect, care should be taken to completely fill the bone defect. One gram of BONITmatrix® is sufficient, to fill a defect of approximately 2 cm². The material is loosely packed into the defect; compression of the material should be avoided as this will reduce the macro porosity.

Other references

BONITmatrix® is gamma-sterilized and is intended for single use. Resterilization is not possible. Products supplied sterile must not be reapplied or applied repeatedly since functionality and sterility can then no longer be guaranteed.

Shelf life

Unopened packets may be stored until the expiry date printed on them. The material must not be used after the expiry date.

Storage conditions

BONITmatrix® must be stored in the sealed original packaging at room temperature. The vial containing the BONITmatrix® must not be taken out of the blister pack until immediately before use.


Time of the publication of the information

January 2021

Other

If you still have questions or comments on the product, you can reach us under the telephone number+49 (0) 381 364669-0 or under info@fmz-gmbh.de.

Legend

	Do not use after the expiry date
	Gamma-sterilized
	Batch number
	For single use only
	Follow the instructions for use
	Manufacturer

CE 0297

GA-BONITmatrix-12012021-0

FMZ GmbH
Charles-Darwin-Ring 3a D-18059 Rostock
Tel. (0381) 364669-0 Fax (0381) 36469810

Freigabe

Datum: 20. JAN. 2021

Name: 